



Universal Base

Product overview USA



Universal Base – engaging

Single unit restoration

35
Uni



Universal Base engaging

Connection	CC			External hex		
	NP	RP	WP	NP	RP	WP
Platform						
1.5 mm collar (article no.)	301101	301102	301103	301107	301108	301109
3 mm collar (article no.)	301104	301105	301106	301110	301111	301112
Post height	4 mm	4 mm	4 mm	4 mm	4 mm	4 mm
Diameter	Ø 4.775	Ø 4.775	Ø 6.515	Ø 4.775	Ø 4.775	Ø 6.515
Screw	37891	37892	37892	29282	29283	29284
Screwdriver	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	37894	37895*	37895*	37894	37895*	31169
Elos Accurate™ Position Locator	IO 2B-A SA	IO 2B-B SA	IO 2B-C SA	IO 6A-A	IO 6A-B	IO 6A-C
Elos Accurate™ Position Locator Kit	IO 2B SA Kit			IO 6A Kit		
Replica for IOS – Elos Accurate™ Analog	PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1	PMA-BRA35-1	PMA-BRA41-1	PMA-BRA51-1
Position Locator DTX Studio	300457	300458	300459	300460	300461	300462
Plaster model implant replica	36697	36698	37879	31158	31159	31160

* 5/pkg

Universal Base Engaging: clinical screw included.

Universal Base – engaging

Single unit restoration

35
Uni



Universal Base engaging

Connection	Tri-Channel			
	NP	RP	WP	6.0
Platform				
3 mm post height (article no.)	301452	301454	301456	301458
5.5 mm post height (article no.)	301453	301455	301457	301459
Abutment collar height	0.5 mm	0.1 mm	0.2 mm	0.4 mm
Diameter	Ø 4.5 mm	Ø 4.5 mm	Ø 5.0 mm	Ø 6.0 mm
Screw	36818	29475	29475	29475
Screwdriver	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	31170	29293*	29293*	29293*
Elos Accurate™ Position Locator	IO 2A-A	IO 2A-B	IO 2A-C	IO 2A-D
Elos Accurate™ Position Locator Kit	IO 2A Kit			
Replica for IOS – Elos Accurate™ Analog	PMA-NBR35-1	PMA-NBR43-1	PMA-NBR50-1	PMA-NBR60-1
Position Locator DTX Studio	300452	300453	300454	300455
Plaster model implant replica	29498	29500	29502	29995

* 5/pkg

Universal Base Engaging: clinical screw included.

Universal Base – non-engaging

Multiple unit restoration – straight design

35
Omni

35
Uni



Universal Base non-engaging – straight design



Connection	CC			External hex			Tri-channel		
	NP	RP	WP	NP	RP	WP	NP	RP	WP
Platform									
Article no.	301362	301363	301364	301367	301368	301369	301370	301371	301372
Abutment collar height	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm
Post height	4 mm	4 mm	4 mm	4 mm	4 mm	4 mm	4 mm	4 mm	4 mm
Diameter	Ø 3.6 mm	Ø 4.2 mm	Ø 5.1 mm	Ø 3.9 mm	Ø 4.4 mm	Ø 5.4 mm	Ø 3.7 mm	Ø 4.5 mm	Ø 5.1 mm
Screw	37367	37606	37606	29282	29283	29284	36818	29475	29475
Screwdriver	Omnigrip	Omnigrip	Omnigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	37374	37607	37607	37894	37895*	31169	31170	29293	29293
Elos Accurate™ Position Locator	IO 2B-A	IO 2B-B	IO 2B-C	IO 6A-A	IO 6A-B	IO 6A-C	IO 2A-A	IO 2A-B	IO 2A-C
Elos Accurate™ Position Locator Kit	IO 2B Kit			IO 6A Kit			IO 2A Kit		
Replica for IOS – Elos Accurate™ Analog	PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1	PMA-BRA35-1	PMA-BRA41-1	PMA-BRA51-1	PMA-NBR35-1	PMA-NBR43-1	PMA-NBR50-1
Position Locator DTX Studio	300470	300471	300472	300463	300464	300465	300466	300467	300468
Plaster model implant replica	36697	36698	37879	31158	31159	31160	29498	29500	29502

* 5/pkg

Universal Base non-engaging: clinical screw included.

Universal Base – multi-unit (non-engaging)

Multiple unit and full arch restoration – straight design

15
O-Mini



Universal Base multi-unit – straight design

Connection	Multi-unit	
	NP / RP ¹	WP ²
Platform	NP / RP ¹	WP ²
Article no.	301365	301366
Abutment collar height	0.3 mm	0.3 mm
Post height	4.5 mm	4.5 mm
Diameter	Ø 5.0 mm	Ø 5.9 mm
Screw	301203	301200
Screwdriver	Omnigrip Mini	Omnigrip Mini
Lab Screw	301201	301202
Elos Accurate™ Position Locator	IO 2C-A	IO 2C-B
Elos Accurate™ Position Locator Kit	IO 2C Kit	
Replica for IOS – Elos Accurate™ Analog	PMA-MUA45-1	PMA-MUA60-1
Position Locator DTX Studio	300473	300474
Plaster model implant replica	38918	31162

¹ Compatible with all Multi-unit Abutments EXCEPT Multi-unit Abutment WP external hex (Brånemark)

² Compatible ONLY with Multi-unit Abutments WP external hex (Brånemark)

Universal Base multi-unit: prosthetic screw included.

On1™ Universal Base

On1™ Universal Base – engaging

Single unit restoration

35
Uni



On1 Universal Abutment – engaging

Connection	On1		
	NP	RP	WP
Platform			
0.3 mm collar (article no.)	301121	301123	301125
1.25 mm collar (article no.)	301122	301124	301126
Post height	3.98 mm	3.98 mm	3.98 mm
Diameter	Ø 4.775 mm	Ø 4.775 mm	Ø 4.775 mm
Screw	38750	38750	38750
Screwdriver	Unigrip	Unigrip	Unigrip
Lab Screw	38751	38751	38751
IOS Healing Cap & scan body (4.5 mm / 5 mm / 6 mm)	38735 / - / 38734	38737 / - / 38736	38739 / 38738 / -
IOS / Desktop analog	38740	38741	38742
Insertion tool	38201	38203	38204
Position Locator DTX Studio	300437	300437	300437
Plaster model implant replica	38731	38732	38733



Delivered with the corresponding On1 Prosthetic Screw.

Screwdrivers & tools

Screwdrivers

Omnigrip™

Omnigrip Screwdriver Machine

20 mm	37379
25 mm	37380
30 mm	37381
35 mm	37382



Omnigrip Screwdriver Manual

20 mm	37376
28 mm	37377
36 mm	37378



Omnigrip™ Mini

Omnigrip Mini Screwdriver Machine

20 mm	300852
28 mm	300853
36 mm	300854



Omnigrip Mini Screwdriver Manual

20 mm	300855
36 mm	300856



Unigrip™

Unigrip Screwdriver Machine

20 mm	29151
25 mm	29152
30 mm	29153
35 mm	29154



Unigrip Screwdriver Manual

20 mm	29148
28 mm	29149
36 mm	29150



Tools

Elos Accurate™ Analog insertion tool and IO driver

Elos Accurate™ Analog Insertion Pin	PMA-AIP-2
Elos Accurate™ Analog Insertion Screw	PMA-AIS-5
Elos Accurate™ Analog Pliers	PMA-AIT-1
Elos Accurate™ IO driver – mini	C18512
Elos Accurate™ IO driver – short	C13485
Elos Accurate™ IO driver – long	C13048

Other tools – Manual Torque Wrench

Manual Torque Wrench Prosthetic	29165
Manual Torque Wrench Adapter	29167



Compatible softwares



For availability of compatible software and download of local production DMEs, please go to nobelbiocare.com/en-us/local-production

ORDER ONLINE

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

store.nobelbiocare.com/us/en/

ORDER BY PHONE

Contact your sales representative, or call our customer service team.

800 322 5001

LIFETIME WARRANTY

On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit:

nobelbiocare.com/warranty



nobelbiocare.com/en-us/universal-base



GMT 83328 US (Rev 00, 03/23) © Nobel Biocare USA, LLC, 2023. All rights reserved. Distributed by: Nobel Biocare.

Disclaimer: Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a dentist or a physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.

Universal Base engaging tri-channel connection Abutments, Universal Base Non-engaging Abutments, Universal Base Multi-Unit Abutments, Elos Accurate™ Position Locators, Elos Accurate™ Position Locator Kits, Replica for IOS Elos Accurate™ Analog and Elos Accurate™ Analog insertion tool and IO driver are manufactured by Elos Medtech Pinol A/S, Engvej 33, 3330 Gørlose, Denmark and distributed by Nobel Biocare.